



# Certificate of Containment Testing

## Containment Testing of Bucket GLBB-bdc-15-100 in a GYROZEN Centrifuge

Report No. 20/110

**Report Prepared For:** GYROZEN Co., Ltd.

**Issue Date:** 12 March 2021

### Test Summary

Bucket GLBB-bdc-15-100 was tested for containment in a GYROZEN centrifuge, using Annex AA of IEC 61010-2-020:2016 (3<sup>rd</sup> Ed.). The sealed buckets were shown to contain a spill.

**Report Written By**

**Name:** Ms Helen Hookway

**Title:** Biosafety Scientist

**Report Authorised By**

**Name:** Mrs Sara Speight

**Title:** Senior Biosafety Scientist



# Certificate of Containment Testing

## Containment Testing of Bucket GLBB-b100-100 in a GYROZEN Centrifuge

Report No. 20/110

**Report Prepared For:** GYROZEN Co., Ltd.

**Issue Date:** 12 March 2021

### Test Summary

Bucket GLBB-b100-100 was tested for containment in a GYROZEN centrifuge, using Annex AA of IEC 61010-2-020:2016 (3<sup>rd</sup> Ed.). The sealed buckets were shown to contain a spill.

#### Report Written By

**Name:** Ms Helen Hookway

**Title:** Biosafety Scientist

#### Report Authorised By

**Name:** Mrs Sara Speight

**Title:** Senior Biosafety Scientist



# Certificate of Containment Testing

## Containment Testing of Rotor GREB-m2.0-24 in a GYROZEN Centrifuge

Report No. 20/040

**Report Prepared For:** GYROZEN Co., Ltd.

**Issue Date:** 27 July 2020

### Test Summary

Rotor GREB-m2.0-24 was containment tested in a GYROZEN centrifuge, using Annex AA of IEC 61010-2-020:2016 (3<sup>rd</sup> Ed.). The sealed rotor was shown to contain a spill.

Report Written By

**Name:** Ms Helen Hookway

**Title:** Biosafety Scientist

Report Authorised By

**Name:** Mrs Sara Speight

**Title:** Senior Biosafety Scientist